2x Realtime PCR Super mix (SYBRgreen,with anti-Taq) 使用说明书

| 产品名称 | 单位 | 货号 |
|--------------------------|--------------|----------|
| 2×Realtime PCR Super mix | 1ml*5 支 | S6013-01 |
| 2×Realtime PCR Super mix | (1ml*5 支)×5 | S6013-05 |
| 2×Realtime PCR Super mix | (1ml*5 支)×10 | S6013-10 |

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2x Realtime PCR Super mix (SYBRgreen, with anti-Taq) User Manual

| Product | Units | Cat.# |
|--------------------------|-----------------|----------|
| 2×Realtime PCR Super mix | 1ml*5tubs | S6013-01 |
| 2×Realtime PCR Super mix | (1ml*5tubes)×5 | S6013-05 |
| 2×Realtime PCR Super mix | (1ml*5tubes)×10 | S6013-10 |

Store at -20°C, protected from light.

Description

S6013 is a Taq DNA polymerase-based 2x Super mix for real-time PCR, which contains all components, except for the primer. This reagent is applicable for intercalation assay with SYBR® Green I. This reagent can be used in glass capillary systems (e.g., LightCycler, Roche Molecular Systems, Inc.). This reagent can be used in a passive reference system (e.g., ABI PRISM® 7700, Applied Biosystems, Inc.). The passive reference dye does not affect any other systems. Hot Start technology with anti-Taq DNA polymerase antibodies enables high specificity and reproducible amplification.

Detection

S6013 can be used in general detection devices, such as: LineGene (Bioer Technology co., Ltd.); S6013 can also be used in detection equipment using glass capillaries or passive reference, such as: LightCycler (Roche Molecular Systems); ABI PRISM® 7000, 7700, and 7900 (Applied Biosystems). Note: The passive reference mode of detectors should be set at "ROX".

Specimen

1. cDNA: Reverse transcription reactions from total or poly (A)+ RNA may be used directly, or after dilution, for real-time PCR. Purified cDNA by phenol/chloroform extraction and ethanol precipitation may also be used. Oligo dT and random primers are suitable for the reverse transcription reaction. Up to 20% of the synthesized cDNA solution from the S6 Super qPCR RT kit

(Code NO. S6012) may be added to the PCR reaction solution directly, without purification.

2. Genomic DNA: Purified DNA, which would be used for general PCR, is also suitable for real-time PCR. In the case of mammalian genomic DNA, 1~10 ng genomic DNA is sufficient for real-time PCR.

Protocol I

Intercalation assay protocol using ABI PRISM® 7700

The following is an intercalator assay protocol to be used with ABI PRISM® 7700. For other detection devices, this protocol may require modification depending on each instruction manual.

1. Preparation of reaction solution

| Component | Volume | Final Concentration |
|--|------------|---------------------|
| PCR grade water | 16 μ l | |
| Realtime PCR Super Mix | 25 μ l | 1x |
| 10pmol/ μ l (10 μ M) Primer #1 | 2 μ l | 0.4 μ M |
| 10pmol/ μ l (10 μ M) Primer #2 | 2 μ l | 0.4 μ M |
| Template DNA | 5 μ l | |
| Total volume | 50 μ l | |

Notes: The primer concentration can be further optimized, if needed. The optimal range for the primers is 0.2~0.6 μ M. In the case of commercially available primers, recommended conditions from those companies should be used.

2. Cycling conditions <3-step cycle> The following condition is recommended:

| | | |
|-------------------|---------------------------------|----------------|
| Pre-denaturation: | 95°C, 30sec. ~1 min. | |
| Denaturation: | 95°C, 15 sec. | |
| Annealing: | 55~65°C, 15 sec. | |
| Extension: | 72°C, 45 sec. (data collection) | (35~40 cycles) |

Melting curve analysis

Notes:

A. The annealing temperature in 3-step cycle should be set to 55~65°C, depending of the primer T_m value.

B. The pre-denaturation condition described above is sufficient for inactivation of the anti-Taq DNA polymerase antibodies used in Hot Start PCR. To prevent unexpected and inappropriate results, do not prolong the pre-denaturation period. Fifteen seconds is also sufficient for denaturation during each cycle.

C. Data collection step should be longer than 30 sec.

Protocol II Intercalation assay protocol using Roche LightCycler™

The following is an intercalator assay protocol to be used with the Roche LightCycler™. In the case of other detection devices, this protocol should be modified accordingly.

A. Preparation of reaction solution

| Component | Volume | Final Concentration |
|--|-------------|---------------------|
| PCR grade water | 6.4 μ l | |
| Realtime PCR Super Mix | 10 μ l | 1x |
| 10 μ mol/ μ l (10 μ M) Primer #1 | 0.8 μ l | 0.4 μ M |
| 10 μ mol/ μ l (10 μ M) Primer #2 | 0.8 μ l | 0.4 μ M |
| Template DNA | 2 μ l | |
| Total volume | 20 μ l | |

Notes: The primer concentration can be further optimized, if needed. The optimal range for primers is 0.2~0.6 μ M. In the case of commercially available primers, recommended conditions from each manual should be followed.

B. Cycling conditions <3-step cycle> The following condition is recommended:

| | | |
|-------------------|---------------------------------|----------------|
| Pre-denaturation: | 95°C, 30 sec. ~1 min. | |
| Denaturation: | 95°C, 5 sec. | |
| Annealing: | 55~65°C, 10 sec. | |
| Extension: | 72°C, 15 sec. (data collection) | (35~40 cycles) |

Melting curve analysis

Notes:

- The annealing temperature can be set to 55~65°C, depending on the primer T_m value.
- The annealing time should be set for 5~20 seconds. Longer annealing time results in increased efficiency, and a shorter time decreases non-specific amplification.
- The pre-denaturation condition described above is sufficient for inactivation of the anti-Taq DNA polymerase antibodies used in Hot Start PCR. To prevent unexpected and inappropriate results, do not prolong the pre-denaturation period. Five seconds is also sufficient for denaturation during each cycle.
- Data collection step should be longer than 10 sec. If commercially available primers or probes are employed, the recommended conditions from each company should be used.

Science Tool